TP de chimie

Chemistry Tutorial

**Ref :
282 070**

Français – p 1

English – p 5

Version : 6010

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1 Description

The ions migration experiment offers the possibility of viewing and identifying certain ions by their colouring characteristics and to monitor their movement, the nature of electrical current and its direction.

Simple to use, rapid and economical in chemical products, the ions migration apparatus is designed for conducting this experiment successfully.

The medium on which the migration takes place is made up of a superposition of sheets of filter paper of the same dimensions as the tank.

The migration of ions is obtained by depositing the compounds containing the coloured ions on the sheets of filter paper soaked in a conductive solution, then connecting the apparatus to a power supply delivering a direct voltage of 6 or 12V.

The approximate time for completing the experiment is 20 to 30 minutes.

Diagram 1

1 – Support box.
2 – Safety socket.
3 – Safety socket.

4 – Removable tank.
5 – Power On/Off indicator.
6 – Removable graphite electrodes.

2 Technical characteristics

Dimensions:

Tank (mm): (L.xl.xh.) 70 × 70 × 10

Box (mm): (L.xh.xp.) 100 × 100 × 50

Electrodes: Removable in graphite
(L.xØ) 100 × 7 mm

Connection to circuit: Banana safety plugs Ø 4 mm.

Power supply: Direct voltage of 6 or 12V.

3 Procedure

3.1 Electrical circuit

The apparatus is set up in a circuit (photo 2) with a power supply delivering a direct voltage of 6 or 12V (type power supply Evolution[®] F6-12). It is possible to mount the apparatus in series with an ammeter for checking the intensity delivered by the power supply.

Photo 2

3.2 Conducting the experiment

- **Migration medium**

The medium on which the migration takes place is made up of a superposition of sheets of filter paper soaked in a conductive solution.

The number of superposed sheets and their level of impregnation determine the intensity of current passing through the migration zone and therefore the time necessary for the migration.

Sheets of filter paper	Filter paper (40×80 mm). Superposition of 3 sheets minimum
Conductive solution	Saturated solution of Na ₂ SO ₄ diluted 3 times ⁽¹⁾ .

It is recommended not to use solutions containing chloride ions. They release chlorine vapours during electrolysis, these vapours are toxic and potentially corrosive for the fastening of electrodes.

The sheets are placed in the location provided on the tank after removing the graphite electrodes. The electrodes are replaced on the clips and thus keep the sheets in place.

The conductive solution is poured on the sheets of filter paper **till they are fully soaked**.

4 Deposit of coloured products

Two different methods are possible:

4.1 Deposit of small crystals.

Draw a line in the middle of the sheet of filter paper (represented in dots), deposit on it the following crystals:

- Copper sulphate (II) Cu²⁺+SO₄²⁻ (pale blue)
- Iron chloride (III) Fe³⁺+ 3Cl⁻ (yellow orange)
- Potassium permanganate K⁺+ MnO₄⁻ (violet)


Diagram 3

4.2 Deposit of coloured solutions.

Deposit a few drops of a mixture of equal volumes of a saturated solution of copper sulphate and a dilute solution of potassium permanganate (the potassium permanganate in solution having a high colouring ability, a very dilute solution may be used so that the deposit does not mask the Cu²⁺ ions).

4.3 Duration of the experiment

After creating the deposits of coloured ions, connect the device to the power supply (the power On/Off indicator is lit).

After about 25 to 30 minutes of operation (depending on the intensity of current in the circuit), we observe while applying the protocol **a**, the migration of coloured ions deposited initially at the centre of the sheet, the anions to the terminal + and the cation to the terminal – of the apparatus (diagram 3).

While depositing the solutions (protocol **b**), we observe the migration of a pale blue (Cu^{2+} ions) over the entire width of the sheet to the terminal – and a violet spot turning to brown (MnO_4^- ions) towards the terminal +.

5 Accessories and maintenance

Sheets of filter paper

The apparatus is delivered with the set of 10 sheets 40×80 mm.

Set of 100 sheets (40×80 mm) Ref. 703 527

Electrodes

Set of 2 replacement graphite electrodes Ref. 283 088

6 Maintenance

Remove and rinse the removable tank and the electrodes with water after each usage.

Clean the fastening clips of the electrodes carefully with an absorbent paper soaked in water and dry them for preventing corrosion.

7 After-sales service

The device is under a 2-year guarantee (except electrodes and the movable tanks), it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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